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1. American Cancer Society, Breast Cancer Resource Center, www.cancer.org, April 2001.
2. China Medical Devices Net, www.zgylqxw.cn/Html/2007-03-10/320313.shtml.
3. Nass SJ, Henderson C, Lashof JC, eds. Mammography and Beyond: Developing Technologies for the Early Detection of Breast Cancer. Washington, DC: National Academy Press. 2001. Prepublication copy:13.
4. Ibid:1.
5. Medical Data International, Inc. U.S. Markets for Diagnostic Oncology Products, 1999-2005. Santa Ana, CA:Medical Data International, Inc. 2000. #RP-481430:2-14.
6. National Academy of Sciences, Institute of Medicine (Recent Reports), www.iom.edu, April 2001.
7. Nass SJ, Henderson C, Lashof JC, eds. Mammography and Beyond: Developing Technologies for the Early Detection of Breast Cancer. Washington, DC: National Academy Press. 2001. Prepublication copy:16.
8. Ibid:14.
9. Ibid:19.
10. Ibid:24.

11. Medical Data International, Inc. U.S. Markets for Diagnostic Oncology Products, 1999-2005. Santa Ana, CA: Medical Data International, Inc. 2000. #RP-481430:2-14.
12. Nass SJ, Henderson C, Lashof JC, eds. *Mammography and Beyond: Developing Technologies for the Early Detection of Breast Cancer*. Washington, DC: National Academy Press. 2001. Prepublication copy:21.
13. Ibid:14.
14. National Cancer Institute, Cancer Information, CancerNet, Types of Cancer, Breast Cancer www.nci.nih.gov, April 2001.
15. Ibid:14.
16. Ibid:1.
17. Folkman J. Tumor angiogenesis: therapeutic implications, *New England Journal of Medicine* 1971; 285:1182-1186.
18. Angiogenesis Foundation, Understanding Angiogenesis, www.angio.org, April 2001.
19. Li WW, Li VW, Tsakayannis D, Casey R, Jaffe M, Atwater LA, eds. *Market Study and Analysis of Angiogenesis- Dependent Diseases*. Cambridge, MA: Angiogenesis Foundation, 2001:17.
20. Ibid:14.
21. Weinberg RA, *One Renegade Cell: How Cancer Begins*. New York, NY: Basic Books. 1998:143-146.
22. Eliceiri BP, Cheresh DA. The role of av integrins during angiogenesis. *Molecular Medicine* 1998;4:741.
23. Angiogenesis Foundation, Understanding Angiogenesis, www.angio.org, April 2001.
24. Li WW, Tumor angiogenesis: molecular pathology, therapeutic targeting and imaging. *Academic Radiology* 2000; 7:800-811.
25. Gasparini G, Brooks PC, Biganzoli E, et al. Vascular integrin avb3: a new prognostic indicator in breast cancer. *Clinical Cancer Research* 1998;4:2625.
26. Feldman F, Habif DV, Fleming RJ, Kanter IE, Seaman WB. Arteriography of the breast. *Radiology* 1967;89:1053-1061.
27. Watt AC, Ackerman LV, Shetty PC, et al. Differentiation between benign and malignant disease of the breast using digital subtraction angiography of the breast. *Cancer* 1985;56:1287-1292.
28. Wells PNT, Halliwell M, Skidmore R, Webb AJ, Woodcock JP. Tumor detection by ultrasound doppler bloodflow signals. *Ultrasound*. 1977;15:231-232.
29. Schoenberger SG, Sutherland CM, Robinson AE. Breast neoplasms: duplex sonographic imaging as an adjunct in diagnosis. *Radiology* 1988;168:665-668.
30. Cosgrove DO, Bamber JC, Davey JB, McKinna JA, Sinnett HD. Color doppler signals from breast tumors. Work in progress. *Radiology* 1990;176:175.
31. Folkman J, Watson K, Ingber D, Hanahan D. Induction of angiogenesis during the transition from hyperplasia to neoplasia. *Nature* 1989;339:58-61.
32. Weidner N, Semple JP, Welch WR, Folkman J. Tumor angiogenesis and metastasis – correlation in invasive breast carcinoma. *New England Journal of Medicine* 1991;324:1-8.
33. Baish JW, Netti PA, Jain RK. Transmural coupling of fluid flow in microcirculatory network and interstitium in tumors. *Microvascular Research* 1997;53:128.
34. Dyachenko A. Dynamic imaging of breast lesions; one dimensional optical model. *Asian Journal of Physics* 2001;10;4:1-18.
35. Boucher Y, Leunig M, Jain RK. Tumor angiogenesis and interstitial hypertension. *Cancer Research* 1996;56:4264.

36. Netti PA, Roberage S, Boucher Y, Baxter LT, Jain RK. Effect of transvascular fluid exchange on pressure—flow relationship in tumors: a proposed mechanism for tumor blood flow heterogeneity. *Microvascular Research* 1996;52:27.
37. Boucher Y, Baxter LT, Jain RK. Interstitial pressure gradients in tissue-isolated and subcutaneous tumors: implications for therapy. *Cancer Research* 1990;50:4478.
38. Boucher Y, Jain RK. Microvascular pressure is the principle driving force for interstitial hypertension in solid tumors. *Cancer Research* 1992;52:5110.
39. Dewhirst MW, Secomb TW, Ong ET, Hsu R, Gross JF. Determination of local oxygen consumption rates in tumors. *Cancer Research* 1994;54:3333.
40. Ertefai S, Profio AE. Spectral transmittance and contrast in breast diaphanography. *Medical Physics* 1985;12:393-400.
41. Profio AE, Navarro GA, Sartorius OW. Scientific basis of breast diaphanography. *Medical Physics* 1989;16:60-65.
42. Eliceiri BP, Cheresh DA. The role of av integrins during angiogenesis. *Molecular Medicine* 1998;4:743.
43. Gasparini G, Brooks PC, Biganzoli E, et al. Vascular integrin av β 3: a new prognostic indicator in breast cancer. *Clinical Cancer Research* 1998;4:2625.
44. Sipkins DA, Cheresh DA, Kazemi MR, Nevin LM, Bednarski MD, Li KCP. Detection of tumor angiogenesis *in vivo* by avb3-targeted magnetic resonance imaging. *Nature Medicine* 1998;4(num 5):623.
45. Cristofanilli,C, et al, *Angiogenesis Modulation in Cancer Research: Novel Clinical Approaches*. Nature Reviews Drug Discovery 2002, June 1(6): 414-26.
46. Sauer G, Deissler H., *Angiogenesis: Prognostic and Therapeutic Implications in Gynecologic and Breast Malignancies*. Current Opinion in Obstetrics and Gynecology; 2003 February, 15(45)-9.
47. Weidner N, Folkman J, Pozza F, Bevilacqua P, Allred EN, Moore DH, Meli S, Gasparini G. *Tumor Angiogenesis: A New Significant and Independent Prognostic Indicator in Early-Stage Breast Carcinoma*. Journal of the Natl. Cancer Institute December 1992.
48. Weidner N, Semple JP, Welch WR, Folkman J., *Tumor Angiogenesis and Metastasis—Correlation in Invasive Breast Carcinoma*. New England Journal of Medicine 1991, Jan. 3, 324(1): 1-8.
49. Gasparini, G., *Clinical Significance of the Determination of Angiogenesis in Human Breast Cancer: Update of the Biological Background and Overview of the Vicenza Studies*. Eur. J. Cancer 1996, 32A: 2485-93.
50. Cutler M. Transillumination as an aid to diagnosis of breast lesions. *Surgery, Gynecology and Obstetrics* 1929;48:721-727.
51. Ertefai S, Profio AE. Spectral transmittance and contrast in breast diaphanography. *Medical Physics* 1985;12:393-400.
52. Profio AE, Navarro GA, Sartorius OW. Scientific basis of breast diaphanography. *Medical Physics* 1989;16:60-65.
53. Peters VG, Wyman DR, Patterson MS, Frank GL. Optical properties of normal and diseased human breast tissues in the visible and near infrared. *Physics in Medicine and Biology* 1990;35:1317-1334.
54. McDonald DM, Choyke PL. *Imaging of Angiogenesis: From Microscope to Clinic*. Nature Medicine 2003, 9(6): 713-725.
55. Ibid.
56. DOBI System data on file, DOBI Medical Systems, LLC, 2001.

57. Alexandra Athanasiou, Daniel Vanel, Corinne Balleyguier, Laure Fournier, Marie Christine Mathieu, Suzette Delaloge and Clasrisse Dromain, Dynamic Optical Breast Imaging: A New Technique to Visualise Breast Vessels: Comparison with Breath MRI and Preliminary Results. European Journal of Radiology 54(2005) 72-79.
58. Kolb TM. Lichy J., Hewhouse J.H., Comparison of the Performance of Screening Mammography, Physical Examination, and Breast US and Evaluation of Factors that Influence Them: An Analysis of 27,825 Patient Evaluations. Radiology 2002, October: 225 (1): 165-75.
59. Kerlikowske K, Carney PA, Geller B, Mandelson MT, Taplin SH, Malvin K, Ernster V, Urban N, Cutter G, Rosenberg R, Ballard-Barbash R. Performance of Screening Mammography Among Women With and Without a First-Degree Relative with Breast Cancer. Annals of Internal Medicine 2000, Dec. 5; 133 (11): 855-63.
60. McDonald DM, Choyke PL. Imaging of Angiogenesis: From Microscope to Clinic. Nature Medicine 2003, 9(6): 713-725.
61. Bartoňková, H., Standara M. , Schneiderová M., THE RESULTS OF DOBI EXAMINATIONS IN MASARYK MEMORIAL CANCER INSTITUTE, Czech Ongological Society CLS JEP, KLINICKA ONKOLOGIE 18 4/2005: 149-151.
62. Gatzemeier W, Scelsi M, Galetti K, Villani L, Tinterri C, Secci A and Costa A, Dynamic Optical Breast Imaing: A non-invasive, adjunctive method to detect breast cancer, San Antonio Breast Cancer Symposium December 2004, Poster 6011.
63. Tabar L, Vitak B, Chen HH et al. Beyond randomized controlled trials : organized mammographic screening substantially reduces breast carcinoma mortality . Cancer 2001;91:1724-1731.
64. Olsen O, Gotzsche PC. Cochrane review on screening for breast cancer with mammography. Lancet 2001;358:1340-1342.
65. Preventive Services Task Force. Screening for breast cancer: recommendations and rationale. Ann Intern Med 2002;137:344-346
66. Kuhl KC, Schrading S, Leutner CC et al. Surveillance of 'high risk' women with proven or suspected familial (hereditary) breast cancer: first mid-term results of a multi-modality clinical screening trial J Clin Oncol 2003;21 (23suppl) :238
67. Smith JA, Andreopoulou E. An overview of the status of imaging screening technology for breast cancer. Ann of Oncol 2004; 15 (Suppl 1):18-26
68. DOBI Technical Paper, DOBI Medical Systems, LLC, August 2001.
69. Laure S. Fournier etc., Dynamic optical breast imaging:: A novel technique to detect and characterize tumor vessels, July 2008.
70. 徐梅等, 动态光学乳腺成像在乳腺病灶鉴别诊断中的价值分析, 人民军医, 43期, 2009年7月.
71. 谭国娟等, 动态光学乳腺成像技术诊断乳腺癌指标初探, 中国医学影像杂志, 34期, 2010年1月.